



# ANTIBACT SUB

## CONCENTRATED DISINFECTANT FOR DIVING EQUIPMENTS

### DESCRIPTION

**ANTIBACT SUB** is a powerful and efficient purifier developed to reduce significantly the proliferation of bacteria in diving material.

**ANTIBACT SUB** decontaminates diving suits, goggles and mouthpieces.

**ANTIBACT SUB** prevents micro-biopic fermentation and destroys nauseous odours.

**ANTIBACT SUB** has a germicide reinforced with a long-life bactericide.

**ANTIBACT SUB** does not damage the skin, eyes and mucous membranes unlike other disinfectants

### TECHNICAL DATA :

Physical form: liquid

pH value (5% Soln) : 7 ( $\pm 1$ )

Density (g/ml): 0.98 ( $\pm 0.02$ )

Solubility : Soluble in water in all proportions

Color : Colourless

Odor : Odorless

### INSTRUCTIONS OF USE :

**ANTIBACT SUB** is used in dilution : 50ml of product for 25 litres of water.

Let diving suit 15 to 30 minutes in the dilution.

No need to rinse.

**Ref. 458M - Bottle 1L**  
**Ref. 358M.5 - Bottle 4L**  
Shelf life : 3 years

Keep out of the reach of children. Reclose the packaging after each use.  
Exposure to sun and extreme temperatures should be avoided. Keep out of direct sunlight.  
Always carry out a test prior to use on the surface to be treated to determine suitability and the appropriate contact time. Before use, read the label safety section for Health and Safety Information or request the safety data sheet to [servicetechnique@mattchem.fr](mailto:servicetechnique@mattchem.fr).

**MATT CHEM<sup>®</sup>**  
**MARINE**

SAS MATT CHEM PRODUCT®  
FRENCH MANUFACTURE  
(92) Bagneux France  
Tel : + 33 (0)1 42 53 73 73  
[www.mattchem.fr](http://www.mattchem.fr) - [info@mattchem.fr](mailto:info@mattchem.fr)  
November 2021